PVC Resuscitation Bag

DE	Beatmungsbeutel aus PVC Gebrauchsanweisung	Seite 3-6
EN	PVC Resuscitation Bag Instructions for use	Page 7-10
CS	PVC resuscitační vak Návod k použití	Strana 11-14
ES	Resucitador manual de PVC Instrucciones de uso	Página 15-18
FR	Insufflateur en PVC Manuel d'utilisation	Page 19-22
IT	Pallone in PVC Istruzioni d'uso	Pagina 23-26
	Symbol	Page 28

Deutsch

ALLGEMEINES

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English

GENERAL

The bag is available in 3 different sizes:

- | | | |
|----------|--------------------|-----------------|
| - Infant | up to 7 kg | Volume: 250 ml |
| - Child | 7–30 kg | Volume: 500 ml |
| - Adult | greater than 30 kg | Volume: 1800 ml |

The resuscitation bag is supplied complete with an intake valve and a non-rebreathing valve (optionally with or without a pressure relief valve).

The resuscitation bag is also available as a set with mask and reservoir bag.

FIELD OF APPLICATION

The resuscitation bag is intended for patients who require resuscitation or temporary breathing support. Ventilation can be performed with air, or with oxygen if required.

SAFETY INFORMATION




- Read and follow the instructions for use carefully before using the product.
- The product may only be used by medical personnel.
- Insufficient ventilation (airflow) can lead to permanent brain damage in patients.
- Before use, the products must be inspected visually for damages (cracks, breakage, etc.) and must be subjected to a functional check (see section "Functional check").
- If ventilation is not possible with this product, another adequate type of ventilation must be initiated.
- It is not permitted to use the resuscitation bag in toxic atmospheres. Furthermore, make sure not to use the resuscitation bag with oxygen supply near smoke or open flame as this may cause fire.
- The product is intended for single use and must not be reprocessed. Reprocessing influences the function of the product. Any reuse entails a potential risk of cross-contamination.
- Avoid contact with oils or lubricants, since this can cause damage to the material.

FUNCTIONAL CHECK

Before use, the products must be inspected visually for damages (cracks, breakage, etc.) and must be subjected to a functional check:

With pressure relief valve

- 1a) Close pressure relief valve. Attach a rebreathing bag (infants and children, 500 ml; adults, 1500 ml) to the non-rebreathing valve. Compress and release the resuscitation bag repeatedly, checking that the rebreathing bag fills. With continued ventilation, the filling and emptying of the rebreathing bag must be observed.


- 1b) Leak test: remove the rebreathing bag. Close the non-rebreathing valve with your thumb and firmly compress the resuscitation bag with closed pressure relief valve. There must be no escape of air. Then open the pressure relief valve and repeat the previous leak test. The pressure relief valve must open and cause the pressure to drop.

Without pressure relief valve

- 2a) Attach a rebreathing bag (infants and children, 500 ml; adults, 1500 ml) to the non-rebreathing valve. Compress and release the resuscitation bag repeatedly, checking that the rebreathing bag fills. With continued ventilation, the filling and emptying of the rebreathing bag must be observed.
- 2b) Leak test: remove the rebreathing bag. Close the non-rebreathing valve with your thumb and firmly compress the resuscitation bag. There must be no escape of air.

With pressure relief valve, adjusted constantly to 40 cmH₂O

- 3a) Attach a rebreathing bag (infants and children, 500 ml; adults, 1500 ml) to the non-rebreathing valve. Compress and release the resuscitation bag repeatedly, checking that the rebreathing bag fills. With continued ventilation, the filling and emptying of the rebreathing bag must be observed.
- 3b) Leak test: remove the rebreathing bag. Close the non-rebreathing valve with your thumb and activate the pressure relief valve with one finger. Compress the resuscitation bag. There must be no escape of air. Then let go of the pressure relief valve and repeat the previous leak test. The pressure relief valve must open and cause the pressure to drop.


If the reservoir bag is used, the following function check must be carried out prior to use:

Attach the reservoir bag and the reservoir valve to the intake valve of the resuscitation bag. Start the flow of air via the oxygen connection into the intake valve. The reservoir bag must fill up. Once the reservoir bag is full, the excess air must flow out via the intake valve. If the resuscitation bag is pumped, the reservoir bag must empty.

Damaged products must not be used.

USE


- Clear mouth and throat of any foreign objects.

- Correctly position the patient so that the airways are open.


- Connect the mask to the non-rebreathing valve.

- Press the mask tightly to the patient's face forming a good seal.

- By compressing and releasing the bag, ventilate the patient. Verify successful ventilation using a suitable method.


- If the resistance is too high during insufflation, check that the airways are free of foreign bodies or correct the head position.
- Watch out for spontaneous respiration!

OXYGEN SUPPLY

Oxygen can be delivered while ventilating (nearly 100% saturation) by attaching a VBM reservoir valve and a VBM reservoir bag to the intake valve of the resuscitation bag and supplying oxygen via the port on the intake valve. If reservoir bags and reservoir valves from other manufacturers are used, the compatibility of the connections must be checked beforehand.

The intake valve has two functions:

1. If the oxygen flow is set too high, any excess oxygen is vented.
2. If the oxygen flow is too low, the missing respiratory gas volume is equalised by taking in air.


DEMAND VALVE

A demand valve can be connected to the intake valve.

PEEP valve

A PEEP valve can be attached to the non-rebreathing valve by means of an adapter.

For further information, refer to the instructions for use of the PEEP valve.

STORAGE AND TRANSPORT CONDITIONS

The maximum shelf life of the product is 5 years.

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international legal regulations.

MATERIAL DATA	
Bag Reservoir bag Oxygen connection tube	PVC (polyvinyl chloride)
Non-rebreathing valve Intake valve	Silicone, PC (polycarbonate)

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