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* For applicability see label on pack

Instructions for the reprocessing of Silkomed products
Disinfectants based on alkylamines and concentrated lyes are unsuitable for use on instruments made of silicone as they cause the material to harden. Suitable disinfectants are those based on aldehydes, active chlorine, active oxygen and quaternary compounds. The exposure times and concentrations recommended by the manufacturer must not be exceeded, otherwise damage may be caused to the material. This silicone product can be sterilized with steam (134 °C / 5 min) or ethylene oxide. For this purpose, the existing filling channels should first be cleared (remove stoppers). Before endotracheal tubes with cuffs made of SILKOMED® are disinfected or cleaned, the air-filling system must be closed in order to prevent liquid penetration. Before endotracheal tubes with cuffs made of SILKOMED® are thermally reprocessed, all air must be removed from the cuff by extraction with a syringe before sterilization takes place. Otherwise there is a risk that the cuff will burst during sterilization.

Caution:
The silicone cuff of the endotracheal tube generally loses air in its inflated condition. For this reason, the cuff volume and pressure should be regularly observed in order to ensure that satisfactory sealing of the trachea is maintained. When endotracheal tubes with cuffs are used during anaesthesia with nitrous oxide, it is vital to check the cuff pressure frequently in order to prevent excessive cuff pressure. When connecting hoses are being thermally recycled, standard connectors must not be sterilized in the connected condition, e.g. with elbow or Y-piece attached. Before reuse, the connection should be checked for firm seating and satisfactory sealing.

Teleflex

 Willy Rüscher GmbH,
Willy-Rüscher-Straße 4–10
71394 Kernen, Germany
Tel. + 49 71 51 40 60
Fax + 49 71 51 40 61 50
www.teleflex.com

93 14 65 Rev. C 08 14 50
en/de/fr/es/it/pt/sv/nl/da/tr/pl
Druckerei Tauber GmbH, Waiblingen
Printed in Germany.
Imprimé en R.F.A.