Deutsch

English

Français

# Gebrauchsanweisung Operating instructions Mode d'emploi

# AtmoSafe

CE

## Inhalt

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Printed by: ATMOS MedizinTechnik GmbH & Co. KG

## Contents

	Page
<b>1.0 Introduction</b>	
1.1 Notes on operating instructions .....	5
1.2 Function .....	7
1.3 Explanation of symbols .....	9
<b>2.0 Safety advice .....</b>	<b>11 - 13</b>
<b>3.0 Intended use.....</b>	<b>14</b>
<b>4.0 Installation and start-up</b>	
4.1 Illustrations.....	15 - 17
<b>5.0 Operation</b>	
5.1 First start-up.....	19
5.2 Installation.....	19
5.2.1 Installation together with a surgical device .....	19
5.2.2 Main filter .....	21
5.2.3 Hose.....	21
5.2.4 Prefilter.....	21
5.2.5 Check on supply voltage.....	21
5.3 Settings .....	23
5.4 Display elements.....	25
5.5 Aspiration .....	25
5.6 Options.....	25
5.7 Service level.....	27 - 29
<b>6.0 Cleaning and maintenance advice</b>	
6.1 Fundamentals of cleaning and disinfection.....	31
6.2 Reprocessing of hoses .....	31
6.3 Cleaning and disinfection of surface of device .....	33
<b>7.0 Maintenance and service .....</b>	<b>35</b>
<b>8.0 Trouble-shooting.....</b>	<b>37</b>
<b>9.0 Accessories and spare parts</b>	
9.1 Accessories.....	39
9.2 Spare parts .....	41
<b>10.0 Technical specifications.....</b>	<b>43</b>
<b>11.0 Disposal.....</b>	<b>45</b>
<b>12.0 Notes on EMC.....</b>	<b>49-51</b>

### General Standard Terms and Conditions

All rights regarding these operating instructions, in particular the right of reproduction and dissemination, and that of translation, are reserved. No part of these operating instructions may be reproduced in any form (by photostat, microfilm, or other processes) without prior written permission of **ATMOS MedizinTechnik GmbH & Co. KG** or by using electronic systems processed, duplicated, or disseminated. The information contained in these operating instructions may be changed or amended without prior announcement and does not represent any liability on the part of **ATMOS MedizinTechnik GmbH & Co. KG**

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## 1.0 Introduction

### 1.1 Notes on operating instructions

- These operating instructions contain important notes on how you can operate the AtmoSafe safely, competently, and effectively. Therefore, it is not only intended for the training of new operating personnel, but also as a reference source. It helps to avoid risks, and also to reduce repair cost and standstill times. Further, it will improve the reliability and life of the device. For these reasons, the operating instructions must always be held available near the unit.

Please read Chapter 2.0 "Safety advice" prior to first start-up.

The following principles must be observed:

**Prudent and careful work is best assurance for protection against accidents!**

Maintenance and repair work may be carried out only by an expert authorised by ATMOS. In case of repairs please insist that only original spare parts are used. Then, you can be assured that operational safety, availability, and value of your device are being preserved.

- The product AtmoSafe has received CE marking in accordance with EU Guideline of the Council on Medical Products 93/42/EWG and meets the basic requirements of Appendix I of this guideline.
- The product AtmoSafe complies with all applicable requirements of the directive 2011/65/EU restricting the use of certain hazardous substances in electrical and electronic equipment ("RoHS").
- The declaration of conformity can be obtained on our website at [www.atmosmed.com](http://www.atmosmed.com).
- The quality management system applied at ATMOS has been certified according to international standards EN ISO 9001 and EN 13485.
- For authorised service work ATMOS is providing service instructions with detailed descriptions of wiring diagrams, setting instructions, and service information.
- Reproduction, also in extracts, only with written permission by ATMOS.

**Abbreviations / symbols used in these operating instructions:**

- Indication of a list
- Subdivision of a list/ activity

The recommended sequence must always be maintained!

☞ Indication of particularly important notes!

☛ Description of the effect of an activity

In the **English** sections, legends are framed like this.

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## 1.2 Function

- The AtmoSafe is an electrically operated medical device, cleaning the air in medically used rooms from fume and gas constituents arising from vaporising of human or animal tissue. This fume arises typically when using laser or electro-surgical devices and is consisting as a rule of water vapour, aerosol, and organic gases.
- The AtmoSafe improves the environmental conditions in the operating theatre:
  - Reduction of the dust load by removing respirable particles
  - Improvement of view
  - Removal of foul-smelling and partly toxic organic gases
  - Evacuation and filtering of hazardous bio-aerosol (filtration of viruses)
- The AtmoSafe draws in the fume-laden air through a device positioned at the surgical application part or through a separate hose, manually held near the application part. By means of a high-efficiency filter the drawn-in air is cleaned of the hazardous matter and returned again into the ambient air.
- Cleaning of the air refers to constituents like aerosol, which is being retained by means of an ULPA high-efficiency *particle filter*, as well as to organic gases, which, amongst other characteristics, feature a smell unpleasant for human beings, retained by a special *gas filter*.

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Fußschalter

### 1.3 Explanation of symbols

Attention, pay attention to operating instructions

Device protection fuse

Equipotential bonding

Device type CF, defi-protected

Alternating voltage

Filter

Delayed stop period

Operation flow

Base flow

Start

Stop

Display: Filter to be changed

Display: Air passage blocked

Unit must not be exposed to explosive anaesthetic gases

Anaesthetic-tested

Protected for use in explosive atmospheres

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## 2.0 Safety advice

- The AtmoSafe is produced in accordance with IEC 601/EN 60601. It is a device of VDE Class of Protection I. It may be connected only to a correctly installed socket outlet with earthing contact.
- Before connecting the device one must check whether the mains voltage specified on the device and mains frequency correspond to the data of the supply network.
- The AtmoSafe may be used only by expert personnel in supervised operation (IEC 601-1/EN 60601-1).
- Following transport at low temperatures the device must be held at ambient temperature for up to four hours before first start-up. When the device has not been acclimatised, it must not be operated.
- The suction hose must never come into direct contact with the evacuation position in order to avoid suction adherence on the tissue.
- Prior to starting the device the power line, accessories, and hoses must be checked for damage. Damaged lines and hoses must be replaced immediately. Prior to use the function of the device is to be tested.
- Disconnection from supply network only by pulling the mains plug! First the plug is to be pulled from the wall socket. Only then the line may be disconnected from the unit.
- Do never touch the plug or the line with wet hands.
- The service environment as specified in the technical data (Chapter 10.0) is to be observed.
- The unit must be positioned on a firm, level, and hard base.
- Operation only in rooms used for medical purposes. The AtmoSafe is not intended for operation within zones subject to explosion hazards (M and G). Areas subject to explosion hazards may arise due to the utilisation of inflammable anaesthetics, skin cleansing agents, and disinfectants. The suction opening of the hose should not lie on the floor. Evacuation of explosive endogenous gases (e.g. methane) from the intestinal tract is also prohibited.
- The foot switch must be suitable for use in areas subject to explosion hazards.  
Foot switches used on operation theatres must have a watertight switching element (IPX 8).
- This product is not re-sterilisable. Repeated reuse of components which are marked with a is forbidden. In case of repeated reuse these components lose their function and there is a high infection risk.
- The AtmoSafe may be operated only in rooms used for medical purposes, but not in areas subject to explosion hazards and in oxygen rich environments.

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## 2.0 Safety advice (Cont.)

- No liquid must penetrate into the device, or be sucked in. When any liquid has penetrated into the device it may be put into operation again only after having been examined by the customer service.
- The AtmoSafe meets the requirements regarding interference immunity of standard **IEC 601-1-2 / EN 60601-1-2** "Electromagnetic compatibility - Medical electric devices".
- ATMOS does not accept liability for personal injury and damage to property when
  - no original ATMOS parts are used,
  - the notes on application in these operating instructions are not observed,
  - installation, re-setting, changes, extensions, and repairs had not been carried out by persons authorised by ATMOS.
- In the surgical invasive sterile area sterile parts only may be used. Hoses, application parts, as well as instruments, must be used either as sterile single-use part or sterilised multiple-use part.
- Power cables must correspond to the applicable national regulations. The suitability for medical applications must be ensured in particular.
- In case of leaks in the system the indication for suction adherence may fail.
- With very low power setting at the HF surgical device auto-activation may possibly not function.
- In case of high electromagnetic and performance-related interference the sensitivity of auto-activation must be reduced.
- These operating instructions correspond to the construction of the device and the status of the applied safety-relevant standards at the time of printing. All rights to change the circuits, processes, names, software programs, and devices given are reserved.
- Please note:  
A medical insulating transformer with earth leakage monitor or any similar safety system acc. to EN 60 601-1 is required, if several devices are connected over one common power supply. The transformer must correspond to the power consumption of all the devices to be connected.

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UK

**Name:** AtmoSafe

**Main functions:** Smoke extraction system for use in conjunction with smoke or gas producing medical appliances, like laser and high-frequency surgical appliances, cauterisers, oscillation saw, and methods like the removal of bone cement in revision endoprosthetics.

**Med. indications/ application:** For extraction and filtering of evaporated burn-off products arising with thermal medical operations by vaporisation of tissue.

For extraction and filtering of aerosols, emitted by oscillation saws (e.g. in autopsy).

For extraction and filtering of vapours arising when two-component adhesives or cement mixtures are mixed and used (e.g. in implant surgical work).

**Specification of the main function:** Aspiration of the air, contaminated with smoke, through a device positioned at the surgical application part or through a separate hose, manually held near the application part. By means of a high-efficiency

filter the drawn-in air is cleaned of the hazardous matter and returned again into the ambient air. Cleaning of the air refers to constituents like aerosols, which is being retained by means of an ULPA high-efficiency particle filter.

**Application organ:** No specific application organ

**Application time:** Short-term use on the patient (< 30 days)

**Application site:** The application site is the clinical, outpatient as well as the practices area. The application of the device may only be performed by medical trained and introduced staff.

**Contraindications:** No application outside of the medical sector. No suction of flammable, corrosive and explosive substances.

**The product is:**  active  not active

**Sterility:** Not necessary

**Single use product / reprocessing:** The device and accessories are partially reusable. For information on reprocessing, cleaning and disinfection see user manual.

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UK

#### 5.1 First start-up

Prior to shipment each AtmoSafe is being inspected by the manufacturer for function and safety. In order to make sure that the device is working safely following transport and installation the following points should be observed: The user should put the device into operation only if the manufacturer or the supplier

1. has carried out a functional test on the device at the place of operation
2. has instructed the person responsible for operating the device on how to handle the device by means of the operating instructions.

*Under any circumstances you must pay attention to the safety advices in Chapter 2.0 prior to first start-up.*

Following transportation at low temperatures (<0°C) the device must be held for up to six hours at ambient temperature before first start-up. When the device has not been acclimatised evacuation cannot be activated. The device is protecting itself.

### 5.2 Installation

The device has to be positioned on a stable base; attention must be paid to the load carrying capacity.

The base must not be soft (foamed material or similar), so that the exhaust air openings are not covered.

#### 5.2.1 Installation together with a surgical device

The power cable of the surgical device must be fitted with an inlet connector for non-heating apparatus according to IEC 320.

This inlet connector for non-heating apparatus is plugged into the auxiliary mains socket (3, Fig. 3) on the rear of the AtmoSafe. The AtmoSafe is now controlling the power input of the surgical device by means of the internal ISA device (Internal Synchronous Activation) and, in this way, effects auto-activation. Auto-activation means that the evacuation mode is switched on as soon as the thermal surgical function of the surgical device (Cut, coag,...) is activated.

Prior to first use of the surgical device auto-activation must be calibrated once. How the AtmoSafe is to be adapted to the surgical device is said in Chapter 5.7.5 Calibration of auto-activation.

Auto-activation of the AtmoSafe works with most surgical devices in the power range 50 W...1400 W (for inst. ERBE ICC 350, ICC 50,...)

Only compatible devices may be connected to the auxiliary mains socket (3, Fig. 3).

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## 5.2.2 Main filter

- After having taken off the packaging and the plugs, please slide the main filter (2, Fig. 1, p. 15) into the filter shaft of the unit until it locks. Take care that the filter cover is on the right side.

When there is no main filter in the unit evacuation cannot be activated.

## 5.2.3 Hose

- After having taken off the packaging the hose ( $\varnothing$  22 mm) is put with one end onto the connecting nipple of the main filter.

## 5.2.4 Prefilter

- Please put the prefilter onto the hose, preferably at the open hose end. The prefilter protects the hose and the main filter against coarse contamination. For more flexibility at the hose end near the field of operation, the filter can be installed between the hose and the main filter. The direction of flow must be taken care of, it is shown on the prefilter. Only a dry and clean filter may be used.

## 5.2.5 Check on supply voltage

- Please check whether the mains voltage data shown on the fuse drawer (2, Fig. 3, p. 17) of the device agree with the values of the mains supply network, and then connect the AtmoSafe to the supply network.
- With surgical applications we recommend to connect the AtmoSafe additionally through the connection with the equipotential bonding of the examination room.

The AtmoSafe is now ready for operation.

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
## 5.1 Settings

### 5.3.1 Presettings general:

The values for operation flow, base flow, and delay can be set individually. By using the scroll button the desired parameter to be changed can be selected. The currently set flow value is shown in the display in percent. The value can be set in % - related to the total rating of the device - by 1%-steps. The currently set delay is shown in the display in seconds. By means of the arrow buttons the value can be increased or reduced. With prolonged pressing of the corresponding arrow button or , setting speed increases, so that also larger changes can be effected quickly.

### 5.3.2 Base flow:

Also during interruptions of work, i.e. when there is no activation signal received from the laser or the HF surgical device or from the foot switch received, the device works with reduced intensity. In order to be able to set the flow intensity of this base flow you should activate the scroll button repeatedly until the indicating light next to the base flow symbol lights up. Now, you press the arrow buttons until the desired value is obtained. This value can be set between 0 and 30 %. The base flow is switched off after a certain period of time (work's setting = 2min.). This time, i.e. the base flow time, can be set on the service level. (See also chapter 5.7 Service level).

### 5.3.3 Operation flow:

During the surgical operation, i.e. when there is an activation signal received from the laser or from the HF surgical device or from the foot switch, the device shall extract the fumes which develop with full intensity. In order to be able to set the flow rate of this operation flow, you should activate the scroll button repeatedly until the indicating light lights up next to the operation flow symbol. Now, you activate the arrow buttons until the desired value is obtained.

### 5.3.4 Delay:

In order to prevent activation and de-activation of operation flow at short intervals during the frequent interruptions between cutting and coagulation phases the flow continues for a certain period of time after the end of the cutting / coagulation procedure (i.e. after the activation signal is off). In order to be able to set the delay period of this operation flow, you should activate the scroll button repeatedly until the indicating light lights up next to the delay symbol. Now, you activate the arrow buttons until the desired value is obtained.

Adjustments from 0...100 sec are possible as well as "0" i.e. unlimited delay. With the adjustment "0" the start button or a connected foot switch has following function (operated once = 'on', operated once again = 'off').

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
## 5.4 Display elements

### 5.4.1 Flow indication

This points out the air resistance arising in the air passage (through the hose and filter system). When all the green bars are lit up, the suction resistance is low. The less bars are lit, the lower the flow. The display permits to draw conclusions, for inst. regarding the degree of blocking of the prefilter.

When the system is completely blocked, all the green bar indicators are off, and the display on the left lights up red.

### 5.4.2 Filter capacity

The display of filter capacity permits quick assessment of the current filter status. The more green bar indicators are lit up, the more capacity there is still left. With full exhaustion of filter capacity the green indicators are off, and the display on the left lights up red.

## 5.5 Evacuation

Please ensure that for invasive treatment the following parts are sterile with each new patient:

- Suction hose incl. suction jets or suction set.
- All parts extending into the operating field, in particular those near or in contact with the patient.
- Switch on the AtmoSafe by means of power switch (❶, Fig.2, p. 16).
- Please ensure that the suction hose is correctly positioned, i.e. that no sensitive parts or swabs or similar, are sucked in.
- Please activate now the flow by pressing the start button or the connected foot switch. The device will continue to work as long as the button is pressed. Thereafter, the flow will continue for the duration of the set delay period (see Chapter 5.3.4).

## 5.6 Options

### 5.6.1 Foot-operated button (Art. No. see accessories)

- Explosion-protected button (AP safety) for switching flow on and off.
  - Connect the foot-operated button.
  - When pressing the foot switch the flow is switched on. After releasing the foot button the device will continue to work for the duration of the set delay period. The function of the foot button is the same as that of the start button.

\* Tip: If the delayed stop (see section 5.3.4) is set in "ON", the foot-operated button has the function of a switch !

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## 5.7 Service level

UK

The service level permits the user to change certain exwork's settings.

You will reach the service level by pressing the SCROLL button (4, Fig. 2, p. 16) and then switching on the device by means of the power switch (1, Fig. 2, p.16). Then in the display will appear.

Basic program - Display "5 0"

Button	Description
↑	Selection of service sub-program (plus)
↓	Selection of service sub-program (minus)
ⓧ	Exit from service mode
◇	Calling the corresponding sub-program

### Sub-programs – overview

Display	Description
5 0	Change brightness of display
5 1	Automatic activation on / off
5 2	Suction adherence indication on / off
5 3	Change base flow period
5 4	Calibration on auto activation (automatic)
5 5	Calibration on auto activation (manual)

### 5.7.1 5 0 Change brightness of display

Display: 0 0 ..... 010 (Brightness value)

Button	Description
↑	Brightness +
↓	Brightness -
ⓧ	Return (incl. store)

*Description:*

Brightness can be adjusted between 0 and 10.

### 5.7.2 5 1 Automatic activation on / off

Display: 1 0 (Activation off) 1 1 (Activation on)

Button	Description
↑	Activation on
↓	Activation off
ⓧ	Return (incl. store)

*Description:*

The activation by means of the integrated current sensor can be switched on and off as required.

### 5.7.3 5 2 Suction adherence indication on / off

Display: 2 0 (Suction adherence indication on)  
2 1 (Suction adherence indication off)

Display	Description
↑	Suction adherence indication on
↓	Suction adherence indication off
ⓧ	Return (incl. store)

*Description:*

Suction adherence indication can be switched on and off as required. With switched off suction adherence indication the device is not switched off and there is also no acoustic warning signal.

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*Description:*

Here the base flow time can be set. The device remains in base flow mode, as determined by this setting. Thereafter the device is switched off. With setting 99 there is no time limit for the base flow.

#### **5.7.5 5.4 Calibration of auto-activation (automatic)**

**Display: 499 .... 400 (Potentiometer setting in %)**

Button	Description
↑	Start of calibration
ⓧ	Cancel

Following auto-calibration manual postprocessing takes place automatically (see 5.7).

*Description:*

Prior to starting auto-calibration the corresponding surgical device must be connected and activated. The set position should be the minimum setting. The device will then search automatically for the corresponding setting value. For validation, automatic calibration is always followed by manual calibration (see 5.7).

#### **5.7.6 5.5 Calibration of auto-activation (manual)**

**Display: 500 .... 599 (Potentiometer setting in %)**

All LEDs (on = activation takes place, off = no activation)

Button	Description
↑	Activation threshold + (less sensitive)
↓	Activation threshold - (more sensitive)
ⓧ	Return (incl. store)

*Description:*

Here sensitivity of automatic activation can be set. For doing this, the corresponding HF device must be connected. The setting must be done in a manner that with activated surgical device the LEDs are lit, and with non-activated surgical device the LEDs are not lit.

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## 6.0 Cleaning and maintenance advice

### 6.1 Fundamentals of cleaning disinfection

- For disinfection the surface and instrument disinfectants listed on page 33 are suitable.

- ☞ There are disinfectants which can cause discoloration of plastic parts, like the filter case, etc.; this however does not effect the function of the parts.
- ☞ The data on concentration and instructions of the manufacturer must always be observed!

### 6.2 Reprocessing of the hoses

- ☞ Please make sure that with each change of the patient the following parts had been disinfected:
- Suction hose incl. suction jet or suction set.
- Washing and disinfecting in an automatic cleaner and disinfector is also possible.  
Thermal disinfection is carried out at 93° C.
  - After disinfecting the parts are to be mounted again (Chapter 5.0 "Operation").

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### 6.3 Cleaning and disinfection of surface of devices

- ☞ It is imperative to pull the power plug before starting to clean and disinfect the surface of the device.
- ☞ Wipe-disinfection:  
Rub the surface of the device with cloth wetted with cleaning or disinfection solution. No liquid must be allowed to penetrate into the device. All the cleaning agents and disinfectants listed below are suitable.
- ☞ When any liquid has penetrated into the device it may be put into operation again only after having been examined by the customer service.
- ☞ Do not use disinfectants containing alcohol.

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### Recommended instrument disinfectants

Disinfectant	Producer
GIGASEPT FF	Schülke & Mayr, Norderstedt
Sekusept PLUS	Henkel, Düsseldorf
Mucozit-T	Merz & Co., Frankfurt/Main

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### Recommended surface disinfectants

Disinfectant	Producer
TERRALIN	Schülke & Mayr, Norderstedt
QUATOHEX	Braun, Melsungen
Incidin Plus	Henkel, Düsseldorf
Pursept-A (Disinfection spray or disinfection cloth)	Merz & Co., Frankfurt/M.

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## 7.0 Maintenance and service

The base device is maintenance-free.

### Filter replacement:

- As described in operating instructions, section 5.2.2 Main filter.

### Preventive maintenance of evacuation system:

- Prior to any use a visual check of the device should be carried out, including the hoses, main and prefilter and mains power cable.  
When using further accessories, like liquid vessel, foot switch, or activation sensors, also these should be subjected to a visual check.
- Replace damaged parts immediately.
- Otherwise, no maintenance at frequent intervals is required.

### Corrective maintenance of system:

- Any changes or repairs may - with due regard to the special requirements for medical products - be carried out only by the manufacturer or by persons expressly authorised by the latter.
- It is mandatory to have a safety-relevant examination carried out in accordance with the respective national regulations (for inst. according to DIN VDE 0751) subsequent to functional problems or a repair.

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## 8.0 Trouble-shooting

AtmoSafe has been subjected to a thorough quality check at the works. If however any problems should occur, you can possibly eliminate these personally, if observing the following notes.

<i>Symptom of problem</i>	<i>Possible causes</i>	<i>Remedy</i>
● Device does not start	<ul style="list-style-type: none"> <li>– Power plug not fitting correctly</li> <li>– Mains voltage not on</li> <li>– Fuse defective</li> </ul>	<ul style="list-style-type: none"> <li>– Check connection at wall socket</li> <li>– Check house fuse</li> <li>– Have fuse replaced by service</li> </ul>
● Insufficient flow	<ul style="list-style-type: none"> <li>– Leaks in hose lines</li> <li>– Main filter is blocked</li> <li>– Prefilter is blocked</li> </ul>	<ul style="list-style-type: none"> <li>– Check hose lines for firm seating</li> <li>– Replace main filter</li> <li>– Replace prefilter</li> </ul>
● Display <b>E 01</b>	<ul style="list-style-type: none"> <li>– Calibration defect of intern pressure sensor</li> </ul>	<ul style="list-style-type: none"> <li>– Service must carry out calibration on site</li> </ul>
● Display <b>E 02</b>	<ul style="list-style-type: none"> <li>– Excess temperature (<math>&gt;69^{\circ}\text{C}</math>) or sub-temperature (<math>&lt;0^{\circ}\text{C}</math>)</li> </ul>	<ul style="list-style-type: none"> <li>– Switch off device and wait until temperature inside the device retuns to normal</li> </ul>
● Display <b>E 03</b>	<ul style="list-style-type: none"> <li>– Data fault on main board</li> </ul>	<ul style="list-style-type: none"> <li>– Switch device off and on again. Should now suction adherence indication no longer work, customer service must carry out calibration on site.</li> </ul>
● Display <b>E 04</b>	<ul style="list-style-type: none"> <li>– Blower does not work</li> </ul>	<ul style="list-style-type: none"> <li>– Switch device off and on again. With renewed defect inform service</li> </ul>

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## 9.0 Accessories and spare parts

### 9.1 Accessories

Connection line for electric equipotential bonding (possibly only in operating theatre).....008.0596.0

⑦ Foot button for AtmoSafe, suitable for use in Zone M (not for operating theatre), explosion-protected, IPX 1 .....445.0061.0

Operation theatre foot button for AtmoSafe, suitable for use in Zone M explosion-protected, IPX 8 .....445.XXXX.X

Standard rail set 25 x 10 mm / 315 mm for lateral mounting to AtmoSafe .....445.0064.0

⑧ Hose holder for insertion into standard rails (25 x 10 mm oder 30 x 10 mm), for air hoses with Ø 22 mm.....445.0066.0

⑨ Surgical handle, with integrated suction channel, with standard international plug-connection for HF surgical devices, incl. air hose Ø 10 mm, 2.5 m, ESU cable length 3 m .....445.0062.0

⑩ Fume evacuation handle of clip-on type, for adaption of mono-polar standard surgical handle, incl. air hose Ø 10 mm and length 2.5 m.....445.0063.0

Air hose Ø 10 mm, of various lengths .....auf Anfrage

Articulated arm with 3 joints for attaching to normal standard rails (25 x 10 mm), extended length approx. 1.3 m, with 5 hose holders for hose with Ø 22 mm.....445.0060.0

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## 9.2 Spare parts

<u>Description</u>	<u>Article No.</u>
Operating instructions.....	445.0000.i
AtmoSafe compl. base package 445.0000.0.....	
Air hose, Ø 22 mm (W), 2,10 m, for single-use, of E.V.A. ....	005.0200.0
Air hose Ø 22 mm, L = 2.70 m made of Hytrel, connecting sockets made of silicone temperature-resistant up to 200 °C .....	005.0201.0
Air hose Ø 22 mm, L = 2.10 m temperature-resistant up to 200 °C .....	005.0203.0
Air hose, internal Ø = 10 mm, L = 1.8 m temperature-resistant up to 200 °C, made of hytrel, connecting sockets made of silicone .....	005.0204.0
① Hose connector straight, Ø 22 mm (M) to Ø 22 mm (M) .....	000.0683.0
② Hose connector straight, Ø 22 mm (M) to Ø 10 mm (M) .....	000.0689.0
③ Hose connector straight, Ø 22 mm (F) to Ø 10 mm (M) .....	000.0688.0
④ Funnel, flat on one side, of PP, with connection in Ø 22 mm (F).....	000.0687.0
⑤ Suction tube, plastic, with ISO cone 22 mm (M) for suction hose Ø 22 mm (F) .....	445.0055.0
⑥ Prefilter, in-line, (HEPA), with ISO connections in Ø 22 mm (M/F).....	445.0044.0
Main filter unit for AtmoSafe, multi-stage gas filter + ULPA particle filter, micro-biocidal-coated compatible with the environment.....	445.0040.0

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## 10.0 Specification

Mains voltage (incl. tolerance)  
arranged by user / installer by inserting the correct fuse drawer  
AtmoSafe WORLD: 100...230VAC ± 10% (total range 90...253VAC) can be  
AtmoSafe EUROPE: 230 V ± 10 %

Mains frequencies (incl. tolerance)	50...60 Hz +/-1%
Mains connection	Device inlet connector for non-heating apparatus
Auxiliary-mains power socket (IEC 320)	120 V / 9.7 A - 230 V / 6.3 A
Power rating	Max. 400 W
Blower air flow (free-flow)	1600 l/min
Device air flow	650 l/min electronically controlled
Activation	Automatic activation, front-panel button, signal at activation input, optionally by foot switch
Operating mode	Suitable for continuous operation
Protective measures (Fuses,...)	The power pack is fitted with fusible cutouts 5 x 20. The motor of the blower is thermally protected.
Fuse (Breaking capacity H)	100 V: T4A (125 V) / 120 V: T4A (125 V) / 230 V: T3,15A (250 V)
Acclimatisation rules	(Prior to start-up): 6 h acclimatisation after transport at low temperature
Maintenance	(Acc. to DIN 31051 Measures for maintaining the desired condition): Filter change and change of single-use articles (by user). The basic device is maintenance-free.
Designed freedom from maintenance	>10 Years (20,000h) (blower, bearings, rubber parts, etc. )
Corrective maintenance	(Acc. to DIN 31051 Measures for restoring the desired condition): Corrective maintenance by ATMOS or by customer service authorised by ATMOS in accordance with the country-specific regulations (in Germany DIN VDE 0751).
Inspection	(Acc. to DIN 31051 Determination of current state): No repeated safety-practical checks as understood by MPG§§22-24. (German law on medical products)
Ambient conditions	(In operation and in transport / storage):
Ambient temperatures	In operation: +10 ... +40°C / in transport and storage: -40 ... +70°C
Air humidity	5...95% (without condensation)
Atmospheric pressure	+700...1060 hPa
Weight:	14 kg
Dimensions (H x W x D )	H 210 mm x W 368 mm x D 410 mm
Noise level	Max. 52 dB(A)@1 m (acc. to ISO 7779)
Interfaces	(Input and output) suction hose connection 7/8"(22 mm), equipotential bonding, activation input, mains connection
Class of protection	I
Applied Part	CF  defibrillator-protected
Main-filter	ULPA, retention = 99.9999 % @ 0.01µm
Canadian Classification	General Hospital & Personal use
Device Group	80FRN
PNC	3
Risk Class	PUMP, INFUSION
Description	Issue Technical Specifications 01.04.2005

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## 11.0 Diposal

### Disposal of base device:

- The packaging consisting of cardboard and foamed polystyrene can be fully recycled or returned to you supplier for repeated use.
- The AtmoSafe does not contain any hazardous goods.
- The material of the case is fully recyclable.
- The components of the product can be disposed like normal electronic scrap. Recyclable materials should be, as far as possible and reasonable, delivered separately to a recycling organisation.

### Disposal of accessories:

- The main filter and other accessories are biologically contaminated material and have to be disposed of by specialist firms in compliance with given regulation.

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## 12.0 Notes on EMC

- Medical electrical equipment is subject to special precautions with regard to EMC and must be installed according to following EMC notes.
- Portable and mobile HF communication facilities can influence medical electrical equipment.
- The use of other accessories, other converters and cables than stated may lead to an increased emission or a reduced interference immunity of the equipment or system.

### 12.1 Guidelines and Manufacturer's Declaration - Emissions

The AtmoSafe is intended for use in the electromagnetic environment specified below. The customer or user of the AtmoSafe should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The AtmoSafe uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The AtmoSafe is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Inapplicable	

The device may not be used directly next to other devices or piled up with other devices. If operation next to or piled with other devices is necessary, please watch the device to check its intended operation in this arrangement.

### 12.2 Guidelines and Manufacturer's Declaration - Immunity

The AtmoSafe is intended for use in the electromagnetic environment specified below. The customer or user of the AtmoSafe should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2	± 6 kV Contact ± 8 kV Air	± 6 kV Contact ± 8 kV Air	Floors should be wood, concrete, or ceramis tile. If floors are synthetic, the relative humidity should be at least 30%.
EFT IEC 61000-4-4	± 2 kV Mains ± 1 kV I/Os	± 2 kV Mains Inapplicable	Mains power quality should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	± 1 kV Differential ± 1 kV Common	± 2 kV Differential ± 1 kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency 50/60 Hz Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Voltage Dips / Dropout IEC 61000-4-11	<p>&lt; 5 % <math>U_T</math> (&gt; 95 % Dip of the <math>U_T</math>) for 0.5 Cycles</p> <p>40 % <math>U_T</math> (60% Dip of the <math>U_T</math>) for 5 Cycles</p> <p>70% <math>U_T</math> (30 % Dip of the <math>U_T</math>) for 25 Cycles</p> <p>&lt; 5 % <math>U_T</math> (&gt;95 % Dip of the <math>U_T</math>) for 5 sec</p>	<p>&lt; 5 % <math>U_T</math> (&gt; 95 % Dip of the <math>U_T</math>) for 0,5 Cycles</p> <p>40 % <math>U_T</math> (60% Dip of the <math>U_T</math>) for 5 Cycles</p> <p>70% <math>U_T</math> (30 % Dip of the <math>U_T</math>) for 25 Cycles</p> <p>&lt; 5 % <math>U_T</math> (&gt;95 % Dip of the <math>U_T</math>) for 5 sec</p>	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AtmoSafe demands continued function even in case of interruptions of the energy supply, it is recommended to supply the AtmoSafe from an uninterruptible current supply or a battery.
NOTE $U_T$ is the mains alternating current prior to application of the test levels.			

### 12.3 Guidelines and Manufacturer's Declaration - Immunity

The AtmoSafe is intended for use in the electromagnetic environment specified below. The customer or user of the AtmoSafe should ensure that it is used in such an environment.

Immunity Test	IEC 60601-Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz	3 V <sub>eff</sub>	Portable and mobile communications equipment should be separated from the AtmoSafe incl. the cables by no less than the distances calculated/listed below.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p><b>Recommended distances:</b></p> $d = (3,5 / V1) * \sqrt{(P)}$ $d = (3,5 / E1) * \sqrt{(P)} \quad 80-800 \text{ MHz}$ $d = (7 / E1) * \sqrt{(P)} \quad 0,8-2,5 \text{ GHz}$ <p>where „P“ is the max. power in watts (W) and D is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed transmitters, as determined by an electromagnetic site (a) survey, should be less than the compliance level (b). Interference may occur in the vicinity of equipment containing the following symbol.</p> 

NOTE 1 With 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines are not like to be applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.

- a The field strength of stationary transmitters, such as base stations of cellular phones and mobile terrain radio equipment, amateur radio transmitters, cbm broadcast and TV stations cannot be predestined exactly. To determine the electromagnetic environment in regard to stationary transmitters, a study of the location is to be considered. If the measured field strength at the location where the AtmoSafe is used exceeds the above compliance level, the AtmoSafe is to be observed to verify the intended use. If abnormal performance characteristics are noted, additional measures might be necessary, e.g. a changed arrangement or another location for the device.
- b Within the frequency range of 150 kHz to 80 MHz the field strength is to be below 3 V/m.

## 12.4 Recommended separations between portable and mobile RF communications equipment and the AtmoSafe

The AtmoSafe is intended for use in electromagnetic environment in which radiated disturbances are controlled. The customer or user of the AtmoSafe can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the AtmoSafe as recommended below, according to the maximum output power of the communications equipment.

Nominal output of the transmitter  W	Separation distance, depending on transmit-frequency m		
	150 kHz to 80 MHz  $d = [ 3.5 / 3 ] \sqrt{P}$	80 MHz to 800 MHz  $d = [ 3.5 / 3 ] \sqrt{P}$	800 MHz to 2,5 GHz  $d = [ 7.0 / 3 ] \sqrt{P}$
0.01	0.12	0.12	0.24
0.1	0.37	0.37	0.74
1	1.2	1.2	2.4
10	3.69	3.69	7.38
100	11.66	11.66	23.32

For transmitters for which the maximum nominal output is not indicated in the above table, the recommended separation distance  $d$  in meters (m) can be determined using the equation belonging to the respective column whereas  $P$  is the maximum nominal output of the transmitter in watts (W) acc. to manufacturer's specification.

NOTE 1 With 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines are not applicable in any case. The propagation of electromagnetic sizes is influenced by absorptions and reflections of buildings, objects and people.

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## 1. General:

Our General Standard Terms and Conditions apply exclusively. Client's terms and conditions which are contrary to or deviate from our General Standard Terms and Conditions are not recognised unless their validity is explicitly confirmed in writing. Our General Standard Terms and Conditions also apply even if we deliver to clients without reservation, in the knowledge of the client's contrary terms and conditions. Our General Standard Terms and Conditions also apply to all future business with that client.

## 2. Proposal - Order Confirmation

Our proposals are subject to change without notice unless otherwise stated in our order confirmation. Each order is only accepted by us following our written order confirmation.

## 3. Orders

Every order requires an exact description of all of our product's details. We assume no liability for errors and damage caused by inaccurate or incomplete ordering details.

## 4. Prices

Unless otherwise stated in the order confirmation, our prices in the order confirmation are ex factory prices and exclude packaging and value added tax. Packaging is charged separately at cost price in the invoice. Value added tax is charged separately in the invoice according to the legal rate on the invoice date. We reserve the right to change prices appropriately should price reductions or increases, especially due to wage settlements, changes in the price of materials or currency fluctuations, be incurred. Proof of such changes will be provided for the client on request.

## 5. Payment Conditions - Balancing

Unless otherwise stated in the order confirmation, our invoices are payable with a 3% discount within 10 days (except for repair and assembly services) or within 21 days from the invoice date net cash; money receipts is decisive for complying with this term. We are entitled to charge interest after the due date at a rate 2% above the relevant basic interest rate of the German Federal Bank. Should the client have payment arrears, we are entitled to charge interest on arrears at a rate 5% above the relevant basic interest rate of the German Federal Bank. Should we be able to prove higher damages due to arrears, we are also entitled to claim these. The client only has the right to balance invoices against its own claims should such claims be confirmed in a court of law or recognised by us. The client does not have the right of retention due to disputed counterclaims.

## 6. Delivery Periods

Fulfilment of our delivery duties requires the punctual and proper fulfilment of the client's duties. The right to defense on the grounds of an unfulfilled contract is reserved. Should the client default in accepting the goods delivery or breach other cooperation duties, we are entitled either to withdraw from the contract or claim compensation for any increased costs incurred up to that time without setting a further deadline. The right to make further claims is reserved. Furthermore, in such cases, the risk of coincidental destruction or a coincidental deterioration in the quality of the delivered goods is transferred to the client in the case of default in accepting such goods or payment arrears. Acts of God or stoppages (due to insufficient supplies of material, industrial disputes etc.) entitle us either to demand an appropriate extension of delivery periods or to partly or entirely dissolve the delivery contract. This does not give the client the right to claim damages. We have fulfilled delivery periods if the delivery goods have left our factory or the client has been informed of the goods' readiness for delivery within such delivery periods. Delivery periods stipulated by the client are not recognised by us unless they form part of our order confirmation. We adhere to legal terms and conditions in cases where, as a result of an undue delay in the delivery for which we are liable, the client is entitled to claim that his interests in a continued fulfilment of the contract have ceased. We also adhere to legal terms and conditions should a delay in delivery be caused by deliberate or grossly negligent action by us or our representatives for which we are responsible. We are also responsible for such actions by our representatives or agents. Should the delivery delay not be caused by our deliberate infringement of contractual duties for which we are responsible, our liability is limited to damage which is regarded as typical for that case. We are liable according to the legal terms and conditions if and in so far as the delivery delay for which we are responsible is caused by an infringement of a substantial contractual duty. In such cases, our liability is also limited to damage which is regarded as typical for that

case. Should the delivery delay be caused by a culpable infringement of non-substantial contractual duties, our client is also entitled to claim a one-off damage compensation worth 3 percentage points of the delivery value of the goods for each week's delay, up to a maximum which is no higher than 15 percentage points of the delivery value of the goods

## 7. Delivery - Familiarisation

In the case of the delivery of devices for the medico-technical industry which require assembly and/or familiarisation for the final customer using specialist trade personnel (such as Ear, Nose and Throat Apparatus and Suction Units), we reserve the right to deliver the goods exclusively to the relevant specialist traders. Should the trader not carry out assembly and/or familiarisation for the final customer, this is carried out by us. In such cases, we reserve the right to charge the client for the additionally created costs. Our specialist traders operate a recording system so that, if necessary, our products can be traced to the final customer. The specialist trader undertakes to immediately report to us all events and risks which must be reported in connection with our products.

## 8. Passage of Risk - Packaging

Unless otherwise stated in our order confirmation, delivery is agreed ex factory. The risk of the goods' damage or loss is therefore transferred to the client as soon as the goods leave the factory or the client is in default of acceptance of the goods. This also applies to cases where we confirm prepaid carriage. Transport packaging and all other packaging according to the packaging regulations is not returnable. Our client is responsible for disposing the packaging at its own cost. Our deliveries are insured by us at the client's expense unless explicitly otherwise agreed. No insurance is arranged in the case of goods which are collected by our clients. In the case of transport damage, claims are only handled by the client receives confirmation of any damage, reduced weight or loss by the shipping company before accepting the delivery.

## 9. Warranty

The client is responsible for examining the delivered goods immediately after receiving them to determine any eventual deficiencies or delivery errors, and to report these immediately. Should the client fulfil this examining and reporting responsibility, and should payment conditions be fulfilled, we shall be liable to the client within the scope of legal regulations. Our period of warranty shall in all cases be two years. Our client can make use of the warranty as follows, so long as he can provide first buyer proof (in the form of an invoice or delivery note) and provided that the product still has the original, unchanged serial number:

- a. We choose whether to fulfil our guarantee by providing repair services free of charge - either on the client's premises or in our factory - or replacing the product. We can also provide these guarantee services through an authorised company;
- b. Should a product be returned to us, the client agrees to send the product in its original or similar packaging, offering the same protection as the original packaging, to our address or any address notified by us.
- c. Our guarantee ceases to apply if changes of any kind have been made to our product, unless such changes have been made by us or a company authorised by us, or have been previously agreed upon in writing by us. Our guarantee also ceases to apply if third parties have carried out repairs to our products or replaced parts thereof. This applies regardless of the fact whether these measures individually or collectively led to a deficiency of the product;
- d. We accept no responsibility for damage defects caused by
  - operational wear and tear;
  - incorrect installation or incorrect or insufficient maintenance;
  - incorrect operation of the product (in contradiction to the operating instructions delivered with the product); - improper use or operating faults; - inappropriate or negligent handling and care, especially with respect to dirt, lime, suction of fluids, inappropriate cleaning and sterilisation;
  - using accessories and/or replacement parts which are not explicitly approved;
  - incorrect assembly and/or initial operation by the client or third parties; - the client's negligence in handling the product; - unacceptable operating conditions, such as humidity, temperatures, the power supply, vibrations.
  - accidents, acts of God, especially lightning, water, fire, public unrest and insufficient ventilation. We are not liable for damage to other objects apart from our product itself, except in the case of any deliberate or grossly negligent actions by us or our representatives or agents. Should no deliberate breach of contract be claimed, our liability

is limited to damage which is regarded as typical for that case. This also applies in the case of our culpable infringement of substantial contractual duties. The indispensable conditions of German Liability Law remain unaffected thereby.

- For second-hand equipment, the period of warranty shall be reduced to a period of twelve months.

## 10. Reservation of Ownership

We retain ownership of our goods until the receipt of all payments arising from the business relationship, including all demands arising from installation orders, subsequent orders, repairs, accessory deliveries and replacement orders. Should we have agreed upon payment on the basis of cheque and bill transactions, the ownership reservation applies until the cheque received by us has been paid in, and does not expire through our credit upon receiving the client's cheque. In the case of a breach of contract by the client, especially payment arrears, we are entitled to repossess our goods. Repossession of our goods represents a withdrawal from the contract, unless explicitly declared in writing by us. We have the right to utilise the product after its repossession, whilst the income from such use is balanced against the client's arrears, after deducting appropriate utilisation costs. The client is responsible for handling the goods with care. Should maintenance and inspection work be necessary, the client must carry these out punctually at his own cost. Our client is entitled to sell the goods he has bought from us in a proper sale transaction. However, he must immediately assign all outstanding claims to the value of the final invoice sum (including value added tax) of our claims to his customers or third parties. The client is entitled to collect this claim even after such assignment. Our right to collect the claim ourselves remains unaffected thereby. We undertake to release the securities to which we are entitled if requested to do so by the client should the realisable value of the our securities be more than 10 percentage points higher than the outstanding claims. We reserve the right to choose the securities to be released.

## 11. Plans and Illustrations

We retain ownership of and copyrights to all plans, illustrations, calculations and other documents which are attached to our proposals. The client must receive explicit written permission before passing these on to third parties. Imitating our legally patented products is forbidden and will be prosecuted.

## 12. Jurisdiction and Place of Performance

Our central office is the place of performance for all disputes in connection with these General Standard Terms and Conditions and the contracts closed with clients under them. This jurisdiction excludes other jurisdiction relating to persons or subject-matter. Furthermore, our client is not entitled to bring charges against us in another court should he file counter-charges, carry out counterbalancing or declare retention. We, however, are entitled to bring charges against our client at their general place of jurisdiction or at another relevant court recognised by German or foreign law. Unless otherwise stated in the order confirmation, our central office is the place of performance.

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