



# BIS™ BRAIN MONITORING SYSTEM: ADVANCED MONITORING FOR IMPROVED OUTCOMES

## POSTOPERATIVE DELIRIUM

In 2015, a multidisciplinary group of scientists and clinicians discovered that a diagnosis of postoperative delirium independently increases the likelihood of institutionalization, 30-day readmission, and extended hospital stay<sup>1</sup>, which may increase cost of care by over 8,300 USD<sup>2</sup>. Delirious patients experience significantly increased odds of dying in the hospital<sup>3,4</sup>, as well as 3 to 6 months after surgery<sup>3,4,5</sup>.

The American Geriatrics Society recently elevated the issue of postoperative delirium; in 2015, they issued a best practice statement concerning risk factors, diagnosis, and perioperative mitigation strategies<sup>6</sup>. The group recommends managing intraoperative exposure to anesthetic agents by employing “processed electroencephalographic monitors of anesthetic depth during intravenous sedation or general anesthesia of older patients to reduce postoperative delirium”<sup>6</sup>, a recommendation echoed shortly after by a consensus group from Brazil<sup>7</sup>.

Mounting evidence associating the use of processed EEG monitors with a reduced risk of postoperative delirium has contributed to this recent interest<sup>5,8,9,10</sup>. Sieber et al. reported over 2.5-times increased odds of postoperative delirium in elderly spinal surgery patients randomized to deep versus light BIS™ monitoring-guided anesthesia<sup>8</sup>. In a larger cohort of elderly general surgery patients, the BIS™ monitoring-guided anesthesia group experienced a 35% relative reduction in the risk of postoperative delirium compared to routine monitoring<sup>9</sup>. A subsequent study in elderly noncardiac surgery patients echoed these results with a 22% relative reduction in the risk of postoperative delirium<sup>5</sup>. More recently, a secondary analysis of 310 cardiac surgery patients found a 33% relative reduction in risk<sup>10</sup>. Taken together in a meta-analysis, these 4 studies indicate reduced odds of developing postoperative delirium by 44% when BIS™ monitoring-guided anesthesia is used<sup>10</sup>. The trials are summarized in the following clinical evidence package.

> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES



## > TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES

## TABLE OF CONTENTS

>	3	SIEBER, ET AL	Sedation depth during spinal anesthesia and the development of postoperative delirium in elderly patients undergoing hip fracture repair.
>	4	CHAN, ET AL	BIS™-guided anesthesia decreases postoperative delirium and cognitive decline.
>	5	RADTKE, ET AL	Monitoring depth of anaesthesia in a randomized trial decreases the rate of postoperative delirium but not postoperative cognitive dysfunction.
>	6	WHITLOCK, ET AL	Postoperative delirium in a substudy of cardiothoracic surgical patients in the BAG-RECALL clinical trial.
>	7	CLINICAL STUDIES	A selection of clinical studies about BIS™ monitoring technology and postoperative delirium.



## Using BIS™ monitoring technology to provide light sedation can help reduce the relative risk of postoperative delirium by more than 50%

### Sedation depth during spinal anesthesia and the development of postoperative delirium in elderly patients undergoing hip fracture repair.

Sieber FE, Zakriya KJ, Gottschalk A, et al.

> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES

<b>Study design</b>	Single-center randomized controlled trial (US)
<b>Arms</b>	Deep sedation: BIS™ monitoring-guided anesthesia titrated to a BIS™ value of approximately 50 Light sedation: BIS™ monitoring-guided anesthesia titrated to a BIS™ value of $\geq 80$
<b>Objective</b>	Compare the risk of postoperative delirium between patients receiving deep and light sedation
<b>N</b>	114
<b>Population</b>	Elderly patients ( $\geq 65$ years) undergoing hip fracture repair with spinal anesthesia
<b>Delirium assessment tool or definition</b>	Confusion Assessment Method
<b>Timing of assessment</b>	Daily in the morning, starting on the second postoperative day until hospital discharge
<b>Results</b>	<ul style="list-style-type: none"><li>▪ 39% of patients developed postoperative delirium</li><li>▪ Deep sedation was associated with<ul style="list-style-type: none"><li>– more propofol (10.2 vs 2.5 mg/kg; <math>P&lt;0.001</math>)</li><li>– less midazolam (1.26 vs 5.53 mg/kg; <math>P=0.04</math>)</li><li>– lower mean BIS™ value (49.9 vs 85.7; <math>P&lt;0.001</math>)</li><li>– longer duration with a BIS™ value <math>&lt; 50</math> (48 vs 4 mins; <math>P&lt;0.001</math>)</li><li>– more than twice the odds of developing postoperative delirium (OR 2.69; <math>P=0.04</math>)</li></ul></li><li>▪ Light sedation was associated with a 52% relative reduction in the risk of postoperative delirium (19 vs 40%; <math>P=0.02</math>)</li><li>▪ The number needed to treat with light sedation to prevent 1 case of postoperative delirium was 4.7 patients</li></ul>
<b>Conclusions</b>	In elderly patients undergoing hip fracture repair under spinal anesthesia, using BIS™ monitoring to titrate anesthesia to lighter levels can help reduce the risk of postoperative delirium by more than 50%



## BIS™-guided anesthetic delivery is associated with 42% lower odds of postoperative delirium

BIS™-guided anesthesia decreases postoperative delirium and cognitive decline.

Chan MT, Cheng BC, Lee TM, Gin T; CODA Trial Group.

> TABLE OF CONTENTS

> SIEBER, ET AL

> **CHAN, ET AL**

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES

<b>Study design</b>	Multicenter randomized controlled trial (China) <i>Cognitive Dysfunction After Anesthesia (CODA) Trial</i>
<b>Arms</b>	Control: Blinded monitoring; routine care to manage anesthetic delivery Intervention: BIS™ monitoring-guided anesthesia titrated to a BIS™ value of 40 – 60
<b>Objective</b>	Determine the association between BIS™ monitoring-guided anesthesia and the risk of postoperative cognitive dysfunction (POCD) and delirium
<b>N</b>	921
<b>Population</b>	Elderly patients (≥ 60 years) undergoing elective major surgery ≥ 2 hours and expected to stay in the hospital ≥ 4 days
<b>Delirium assessment tool or definition</b>	Acute fluctuating course of inattention, AND disorganized thinking OR altered level of consciousness
<b>Timing of assessment</b>	Daily in the morning, starting on the second postoperative day until hospital discharge
<b>Results</b>	<ul style="list-style-type: none"> <li>▪ 20% of patients developed postoperative delirium</li> <li>▪ BIS™ monitoring-guided anesthesia was associated with <ul style="list-style-type: none"> <li>– less propofol and volatile anesthetic gas use</li> <li>– higher mean BIS™ values (53.2 vs 38.6; P&lt;0.001)</li> <li>– shorter duration with a BIS™ value &lt; 40 (7.2 vs 22.8 mins; P&lt;0.001)</li> </ul> </li> <li>▪ Fewer patients experienced postoperative cognitive issues in the BIS™ monitoring group <ul style="list-style-type: none"> <li>– 38% lower odds of developing POCD at 3 months (P=0.02)</li> <li>– 35% relative reduction in the risk of postoperative delirium (15.6 vs 24.1%; P=0.01)</li> <li>– 42% lower odds of developing postoperative delirium (P=0.01)</li> </ul> </li> <li>▪ In 1000 elderly patients undergoing major surgery, use of BIS™ monitoring-guided anesthesia is expected to prevent 23 cases of POCD and 83 patients of postoperative delirium</li> </ul>
<b>Conclusions</b>	Use of BIS™ monitoring-guided anesthesia reduced the time spent with BIS™ values < 40, subsequently reducing the risk of postoperative cognitive issues



## Use of BIS™ monitoring technology is associated with a reduced incidence of deep anesthesia and 22% lower relative risk of postoperative delirium

Monitoring depth of anaesthesia in a randomized trial decreases the rate of postoperative delirium but not postoperative cognitive dysfunction.

Radtko FM, Franck M, Lendner J, Krüger S, Wernecke KD, Spies CD.

> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> **RADTKE, ET AL**

> WHITLOCK, ET AL

> CLINICAL STUDIES

<b>Study design</b>	Single-center randomized controlled trial (Germany)
<b>Arms</b>	Control: Blinded monitoring; routine care to manage anesthetic delivery Intervention: BIS™ monitoring-guided anesthesia
<b>Objective</b>	Determine the association between BIS™ monitoring-guided anesthesia and the risk of postoperative delirium
<b>N</b>	1,155
<b>Population</b>	Elderly patients (≥ 60 years) undergoing elective major surgery ≥ 1 hour
<b>Delirium assessment tool or definition</b>	Diagnostic and Statistical Manual of Mental Disorders (DSM IV)
<b>Timing of assessment</b>	Twice daily starting on the first postoperative day through the seventh day
<b>Results</b>	<ul style="list-style-type: none"> <li>▪ 18.8% of patients developed postoperative delirium</li> <li>▪ BIS™ monitoring-guided anesthesia was associated with a <ul style="list-style-type: none"> <li>– lower number of average BIS™ values &lt; 20 (P=0.04)</li> <li>– 22% relative reduction in postoperative delirium (16.7 vs 24.1%; P=0.036)</li> <li>– trend towards lower risk of POCD at 7 days (18.1 vs 23.9%; P=0.062)</li> </ul> </li> <li>▪ The percentage of BIS™ values &lt; 20 was independently associated with a higher odds of postoperative delirium (1.027; P=0.006)</li> </ul>
<b>Conclusions</b>	BIS™ monitoring was associated with a reduced incidence of low BIS™ values and a reduced risk of postoperative delirium





## BIS™ monitoring-guided anesthesia is associated with a 44% reduction in the odds of postoperative delirium across 4 randomized controlled trials

### Postoperative delirium in a substudy of cardiothoracic surgical patients in the BAG-RECALL clinical trial.

Whitlock EL, Torres BA, Lin N, et al.

> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES

<b>Study design</b>	Single-center (US) randomized controlled trial Secondary analysis of 6,100- patient multicenter international <i>BIS or Anesthetic Gas to Reduce Explicit Recall</i> (BAG-RECALL) Trial
<b>Arms</b>	Control: Blinded BIS™ monitoring; end-tidal anesthetic concentration (ETAC) monitoring to manage anesthetic delivery to an age-adjusted minimum alveolar concentration of 0.7 – 1.3 Intervention: BIS™ monitoring-guided anesthesia titrated to a BIS™ value of 40 – 60
<b>Objective</b>	Determine the association between BIS™ monitoring-guided anesthesia and intraoperative awareness (the postoperative delirium endpoint was part of a prespecified single-center substudy)
<b>N</b>	310
<b>Population</b>	Cardiac and/or thoracic surgery patients at high risk for intraoperative awareness
<b>Delirium assessment tool or definition</b>	Confusion Assessment Method for the Intensive Care Unit
<b>Timing of assessment</b>	Twice daily until the tenth postoperative day or until ICU discharge
<b>Results</b>	<ul style="list-style-type: none"><li>▪ 23.5% of patients developed postoperative delirium</li><li>▪ BIS™ monitoring-guided anesthesia was associated with a trend towards a 33% relative reduced risk of postoperative delirium (18.8 vs 28.0%; P=0.058)</li><li>▪ When taken together, a <i>posthoc</i> meta-analysis combining the results of this study with three others<sup>5,8,9</sup> found that the use of BIS™ monitoring-guided anesthesia was associated with significant 44% reduction in the odds of postoperative delirium</li></ul>
<b>Conclusions</b>	Although this study did not show a significant association between BIS™ monitoring and a reduced risk of postoperative delirium compared to ETAC, a <i>posthoc</i> meta-analysis including these results showed a significant reduction in the odds of developing postoperative delirium.



## A selection of clinical studies about BIS™ monitoring technology and postoperative delirium

Fritz BA, Kalarickal PL, Maybrier HR, Muench MR, Dearth D, Chen Y, Escallier KE, Ben Abdallah A, Lin N, Avidan MS. Intraoperative electroencephalogram suppression predicts postoperative delirium. *Anesth Analg*. 2015 [Epub ahead of print]

Soehle M, Dittmann A, Ellerkmann RK, Baumgarten G, Putensen C, Guenther U. Intraoperative burst suppression is associated with postoperative delirium following cardiac surgery: a prospective, observational study. *BMC Anesthesiol*. 2015;15:61.

Seo JS, Park SW, Lee YS, Chung C, Kim YB. Risk factors for delirium after spine surgery in elderly patients. *J Korean Neurosurg Soc*. 2014;56(1):28-33.

Brown CH 4th, Azman AS, Gottschalk A, Mears SC, Sieber FE. Sedation depth during spinal anesthesia and survival in elderly patients undergoing hip fracture repair. *Anesth Analg*. 2014;118(5):977-80.

Berger M, Nadler J, Mathew JP. Preventing delirium after cardiothoracic surgery: provocative but preliminary evidence for bispectral index monitoring. *Anesth Analg*. 2014;118(4):706-7.

Whitlock EL, Torres BA, Lin N, Helsten DL, Nadelson MR, Mashour GA, Avidan MS. Postoperative delirium in a substudy of cardiothoracic surgical patients in the BAG-RECALL clinical trial. *Anesth Analg*. 2014;118(4):809-17.

Radtke FM, Franck M, Lendner J, Krüger S, Wernecke KD, Spies CD. Monitoring depth of anaesthesia in a randomized trial decreases the rate of postoperative delirium but not postoperative cognitive dysfunction. *Br J Anaesth*. 2013;110 Suppl 1:i98-105.

Chan MT, Cheng BC, Lee TM, Gin T; CODA Trial Group. BIS-guided anesthesia decreases postoperative delirium and cognitive decline. *J Neurosurg Anesthesiol*. 2013;25(1):33-42.

Santarpino G, Fasol R, Sirch J, Ackermann B, Pfeiffer S, Fischlein T. Impact of bispectral index monitoring on postoperative delirium in patients undergoing aortic surgery. *HSR Proc Intensive Care Cardiovasc Anesth*. 2011;3(1):47-58.

Plaschke K, Fichtenkamm P, Schramm C, Hauth S, Martin E, Verch M, Karck M, Kopitz J. Early postoperative delirium after open-heart cardiac surgery is associated with decreased bispectral EEG and increased cortisol and interleukin-6. *Intensive Care Med*. 2010;36(12):2081-9.

Sieber FE, Zakriya KJ, Gottschalk A, Blute MR, Lee HB, Rosenberg PB, Mears SC. Sedation depth during spinal anesthesia and the development of postoperative delirium in elderly patients undergoing hip fracture repair. *Mayo Clin Proc*. 2010;85(1):18-26.

> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES



> TABLE OF CONTENTS

> SIEBER, ET AL

> CHAN, ET AL

> RADTKE, ET AL

> WHITLOCK, ET AL

> CLINICAL STUDIES

1. Gleason LJ, Schmitt EM, Kosar CM, Tabloski P, Saczynski JS, Robinson T, Cooper Z, Rogers SO Jr, Jones RN, Marcantonio ER, Inouye SK. Effect of delirium and other major complications on outcomes after elective surgery in older adults. *JAMA Surg.* 2015; 9:1-7.
2. Zywielski MG, Hurley RT, Perruccio AV, Hancock-Howard RL, Coyte PC, Rampersaud YR. Health economic implications of perioperative delirium in older patients after surgery for a fragility hip fracture. *J Bone Joint Surg Am.* 2015;97(10):829-36.
3. Abelha FJ, Luis C, Veiga D, Parente D, Fernandes V, Santos P, Botelho M, Santos A, Santos C. Outcome and quality of life in patients with postoperative delirium during an ICU stay following major surgery. *Crit Care.* 2013;17(5):R257.
4. Veiga D, Luis C, Parente D, Fernandes V, Botelho M, Santos P, Abelha F. Postoperative delirium in intensive care patients: risk factors and outcome. *Rev Bras Anesthesiol.* 2012;62(4):469-83.
5. Radtke FM, Franck M, Lendner J, Krüger S, Wernecke KD, Spies CD. Monitoring depth of anaesthesia in a randomized trial decreases the rate of postoperative delirium but not postoperative cognitive dysfunction. *Br J Anaesth.* 2013;110 Suppl 1:i98-105.
6. American Geriatrics Society Expert Panel on Postoperative Delirium in Older Adults. Postoperative delirium in older adults: best practice statement from the American Geriatrics Society. *J Am Coll Surg.* 2015;220(2):136-48.e1.
7. Nunes RR, Fonseca NM, Simões CM, Rosa DM, Silva ED, Cavalcante SL, Lopes CG, Stefani LC. Brazilian consensus on anesthetic depth monitoring. *Braz J Anesthesiol.* 2015;65(6):427-36.
8. Sieber FE, Zakriya KJ, Gottschalk A, Blute MR, Lee HB, Rosenberg PB, Mears SC. Sedation depth during spinal anesthesia and the development of postoperative delirium in elderly patients undergoing hip fracture repair. *Mayo Clin Proc.* 2010;85(1):18-26.
9. Chan MT, Cheng BC, Lee TM, Gin T; CODA Trial Group. BIS-guided anesthesia decreases postoperative delirium and cognitive decline. *J Neurosurg Anesthesiol.* 2013;25(1):33-42.
10. Whitlock EL, Torres BA, Lin N, Helsten DL, Nadelson MR, Mashour GA, Avidan MS. Postoperative delirium in a substudy of cardiothoracic surgical patients in the BAG-RECALL clinical trial. *Anesth Analg.* 2014;118(4):809-17.

# Medtronic



Use scan app to read

For more information, please visit  
[medtronic.eu/product-catalog](http://medtronic.eu/product-catalog)

**IMPORTANT:** Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 16-eu-bis-clinical-evidence-interactive-guide-792054