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BeneHeart Trainer

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Direction for Use

Product Introduction

The BeneHeart AED (Automated External Defibrillator) Trainer is designed specific for training first-aid personnel before they use BeneHeart D1 AED. The trainer fully simulates the BeneHeart D1 defibrillator in the training guide and drill program with a view to help first-aid personnel to get familiar with the AED and grasp required basic operation skills when using it in the first-aid process.

Note: Before using an AED clinically, first-aid personnel must first complete required first-aid training courses so that they can cope with various situations calmly in case of an accident. Such first-aid trainings also include studying the use of the AED trainer. Well-trained first-aid personnel must know very well how to assess the patient's conditions such as loss of consciousness, respiratory arrest, cardiac arrest etc. They are also proficient in CPR operation procedures and related medical procedures.

Product Function Features

BeneHeart trainer (AED trainer) can simulate shock but cannot output high energy, thus ensuring its safety during use. (Note: high-voltage shock is applied in real AED application, be sure to keep this in mind for operation safety).

The BeneHeart Trainer (AED trainer) has preset 6 training scenarios (refer to the detailed description of simulated cases on Page 13 in the manual), which are demonstrated by pictures and sounds; each scenario simulates sudden cardiac arrest of a real heart, which will help its users master the AED operation skills. These scenarios conform to the training standards required by AHA CPR guidelines. Please refer to the local medical guide and emergency medicine for other requirements of the training.


Direction for Use

Product Composition and Button Description

The product consists of main unit, remote control, electrode wire, pads, operator's manual and carry case.

Main Unit

1.1 Front View of the Main Unit

- 1. Pads connector
It is used to connect the pads jack.
- 2. Status indicator
Green: Power indicator.
- 3. Display screen
- 4. Soft key area
Correspond to hotkeys on the screen. Under different interfaces, the functions of the same soft key are different.
- 5. Shock button
Press this button to deliver a shock.
- 6. Handle

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Direction for Use


7. On/Off button

Press this button to power on or off the trainer.

8. Remote control receiver tube

When operating with a remote control, point it to the receiver tube.

1.2 Side View


1. Loud speaker

2. USB interface

Direction for Use

1.3 Rear View

1. Pads compartment

2. Battery compartment

3. Release button for the battery compartment

Press down this button and slide the battery cover to the right to open the battery compartment.

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

Direction for Use

Screen Display

1. Training scenario and simulation mode

There are a total of six training scenarios and eight simulation modes.

2. Number of shocks delivered

3. Timing area of the system

This area shows the operating time of the AED trainer after it is turned on.

4. Patient type (Adult/ Pediatric)

5. Volume icon

6. Battery status indicator

It indicates the battery status.

7. Treatment information area

8. Soft key operation area

Its soft key names correspond to the following soft buttons, which will change with the work mode and the interface view, when the soft button name is blank, it indicates the soft button is inactive currently.

Soft Button Icons


The following icons may appear at the soft button operation area, which meanings are as follows:

- Return to the previous page
- Adjust the parameter by +
- Adjust the parameter by -
- Move the option cursor
- Switching between the Adult/Pediatric Mode
- Power off

Direction for Use

- Setting
- Switch the breathing- cardiac compression rate
- Audio language soft key
Switch the language of audio prompts. The symbol changes if the system language is changed. The larger symbol indicates the current language while the smaller one indicates the target language.


Remote Control

- 1. Scenario selection: Select from 1 to 6.
- 2. Play/pause button: When the trainer is playing, press this button, the trainer enters the pause state; when the trainer is under the pause state, press the button, it will continue its play.
Note: press and hold the Play/pause button can switch the language. When the trainer is playing, press and hold the Play/pause button for about 2 seconds till the language icon is switched, the language will be changed successfully.
- 3. Simulate the scenario that the pads have been connected properly
Press the button during scenario training process, the upper left corner of the main unit screen will display “OK”, it enters the scenario of simulating good pads contact.

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Direction for Use

4. Simulate “shock advised” scenario (shockable rhythm)

Press the button during scenario training process, the upper left corner of the main unit screen will display “⚡”, it enters the scenario of simulating shock.

5. Simulate the scenario of low battery

Press the button during scenario training process, the upper left corner of the main unit screen will display “🔋”, it enters the scenario of simulating low battery.

6. Simulate the scenario of motion artifact

Press the button during scenario training process, the upper left corner of the main unit screen will display “👤”, it enters the scenario of “motion artifact”.

7. Simulate the scenario that the pads are not connected properly

Press the button during scenario training process, the upper left corner of the main unit screen will display “✖️”, it enters the scenario of simulating poor contact of pads.

8. Simulate “no shock advised” scenario

Press the button during scenario training process, the upper left corner of the main unit screen will display “🚫”, it enters the scenario of simulating no-shock.

9. Simulate the scenario that the battery power runs out

Press the button during scenario training process, the upper left corner of the main unit screen will display “🔌”, it enters the scenario of simulating the battery power is running out.


10. Simulate the scenario of unit error

Press the button during scenario training process, the upper left corner of the main unit screen will display “🔑”, it enters the scenario of simulating the unit error.


11. Volume +/-: High/Medium/Low/Mute adjustable.

Other Components

Pads & electrode wires

Operator's manual

Carry case

Direction for Use

Basic Operations and Settings

Battery installation

- 1) Press the battery compartment release button and push the battery compartment cover to the right as per the direction pointed by the arrow shown below.

- 2) Install 8 AA batteries properly, pay attention to the battery polarity. (Note: Batteries used by a true clinical AED are disposable manganese dioxide lithium batteries)
- 3) Close the battery cover of the battery compartment and push it leftward until the battery cover gives out an sound showing it is in position.

Connecting Pads to Electrode Wires

1. Connect the pads and the electrode wires properly. The adult electrode pad must be connected to the adult electrode wire (with orange connector) and the pediatric electrode pad must be connected to the pediatric electrode wire (with green connector). (Note: the electrode wires and pads of the BeneHeart D1 AED are connected together, which do not need additional connection)
2. Plug the electrode wire into the pads jack on the upper left corner of the machine.
3. Place the pads into the pads compartment properly at the back of the unit.

Training Pads

Pads for training are similar to those used in a clinical AED unit, however, their external dimensions and materials are different and cannot conduct current. After the electrode wires are connected to pads, it may be used with any CPR manikin; when inserting the electrode wire terminal into the trainer jack, the trainer will demonstrate the process according to the training procedure.

Note: The training pads are for training use only and must be not used with the real clinical AED.

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


Pediatric Electrode Pad

AED Trainer Setting

- 1. Press the Power button to turn on the trainer.
- 2. Press the Power button again to enter the power-off interface and the parameter-setting interface (see the figure below).


- 2.1 Press the button to turn off the trainer.
- 2.2 Press the button to return to the current scenario.
- 2.3 Press the button to set parameters.

Direction for Use

Parameter setting includes: volume, scene, CPR pattern, charge level and language. Press the >> button to choose parameters such as scene, and then press the Up/Down button to choose different scenes such as 1-6.


Special Tips: this AED Trainer comes with multi-languages function. The default languages are Chinese and English. Users can choose one more language as the third language, and then they can change language among these three languages during operation. For example: users can press the >> button to choose parameter 地球 “Language setting”, and then press the Up/Down button to choose different language such as Russian, after that, they can change language among Chinese-English-Russian during operation. If the third language which the user chose is Chinese or English, then only Chinese and English can be changed during operation.

Note: Set parameters will be saved automatically rather than by pressing any key. Users may carry out other operations such as turn off the trainer or return to the current scenario after setting parameters properly. The trainer will default the parameters set last time when it is powered on next time; For example, if Scene 6 was set last time, the trainer will enter Scenario 6 automatically when it is powered on.

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
Product Operating Procedures


1. Turn on the trainer

Press the Power button to turn on the trainer. The screen will show the icon “call emergency” (as show below) for 2 seconds. Then the AED Trainer will detect the pads connection. If the pads have not been connected and plugged into the pads connector properly, the trainer will keep on prompting its user to connect the pads and the electrode wire plugs properly; if the adult pads are connected, the AED trainer will enter the Adult Mode; if the pediatric pads are connected, the AED trainer will enter the Pediatric Mode.

2. Apply pads

Take out the pads for training and apply them to a CPR manikin body according to the following positions.

3. Carry out operations according to sound and prompt information.

- ◆ If a shockable rhythm is detected, the AED trainer will prompt the user to press the Shock button to simulate a shock.
- ◆ If no shock is advised, the AED Trainer will enter the CPR pattern or keep on analyzing according to the current configurations

4. Select training scenario

The AED trainer has set 6 basic scenarios and 8 optional simulation scenes. Select different scenarios with its remote control for training.

Direction for Use

Product Storage Method

Pads Maintenance Method

- Pads must be kept clean to avoid dust, oily substances, sticky articles and other contaminants, otherwise, pads viscosity will decline.
- Stick the pads to their original covers rather than to other articles randomly for preservation to avoid damaging them.
- When pads viscosity declines due to contamination, clean directly with a small amount of water slightly, which viscosity will be recovered temporarily after being fully dried. If too much water is used for cleaning, pads viscosity will decline.
- Do not wipe the sticky surface with tissue or cloth.
- Do not scrape the sticky surface with nails or brushes.
- Do not wash pads frequently and do not clean pads with cleaning agent or hot water.
- Keep the pads away from direct sunlight, high temperature, moisture, dust and corrosive gas.

Main Unit Maintenance Method

The main unit must be cleaned on a regular basis. At places with heavy environmental pollution or strong sandstorm, the main unit must be cleaned more frequently. Optional cleaning agents include hydrogen peroxide solution (3%) and ethanol (70%).

When cleaning the trainer:

1. Turn off the power and remove batteries.
2. Cleanse the display screen with soft cotton balls soaked with proper amount of cleaning agent.
3. Cleanse the trainer surface with a soft cloth soaked with proper amount of cleaning agent.
4. Remove excess cleaning agent with a dry cloth when necessary.
5. Place the trainer in a ventilated cool environment to let it dry.

Considerations in Trainer Keeping:

- Do not immerse the trainer into liquid.
- Do not pour any liquids onto the trainer or its accessories.
- Do not allow liquid to enter the enclosure.
- Do not use any abrasive material such as steel wool or polishing compound or any strong solvent such as acetone or cleaning agent containing acetone component for cleaning and disinfection.
- Keep the trainer away from direct sunlight, high temperature, moisture, dust or corrosive gas.
- Keep the main unit, pads, electrode wires, remote control and manual into the carry case.

※ Be sure to obey the above considerations to use the trainer correctly.

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

Direction for Use

6 kinds of Scenarios Preset

| Simulated Scenarios | Description | Operation |
|---------------------|--|--|
| 1 | Shockable rhythm— Defibrillation succeeds in one shock | Shockable rhythm One shock Normal rhythm |
| 2 | Shockable rhythm— Defibrillation succeeds in two shocks | Shockable rhythm Two shocks Normal rhythm |
| 3 | Shockable rhythm— Defibrillation succeeds in three shocks | Shockable rhythm Three shocks Normal rhythm |
| 4 | Shockable rhythm – Defibrillation succeeds in one shock Shockable rhythm Defibrillation succeeds in one shock | Shockable rhythm One shock Normal rhythm Shockable rhythm One shock Normal rhythm |
| 5 | Shockable rhythm | Shockable rhythm throughout the whole course |
| 6 | Normal rhythm | Normal rhythm throughout the whole course |

Direction for Use

Product Specifications and Configurations

Product model: BeneHeart Trainer
Product name: AED Trainer
Power supply: DC12V, 8 AA batteries
External dimension: 288X210X80mm
Product weight: 1.4±0.2kg (Standard Edition)
2.8±0.2kg (Weight Increase Edition)
Maximum working current: ≤500mA
Off current: ≤10μA
Note: Support chargeable AA batteries, rechargeable batteries and charger must be purchased by the customer.

Remote control specifications
Power supply: DC 3V, 2 AAA batteries
Static working current: < 5μA
Maximum working current: < 10mA

Product Configurations:

| SN. | Contents | Qty | Unit | Remarks |
|-----|-------------------------|-----|------|--|
| 1 | Main unit | 1 | PC | |
| 2 | Remote control | 1 | PC | |
| 3 | Electrode wire | 2 | PC | 1pc adult electrode wire 1pc pediatric electrode wire |
| 4 | Adult electrode pad | 2 | Pair | Can be used repeatedly |
| 5 | Pediatric electrode pad | 1 | Pair | Can be used repeatedly |
| 6 | Operator's manual | 1 | PC | |
| 7 | Carry case | 1 | PC | |

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