# 551B Veterinary SpO<sub>2</sub> Sensor

## Instructions for Use

### Statement

**SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.** (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden. This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures animal and operator safety.

Mindray is responsible for safety, reliability and performance of this product only in the condition that:

- This product is operated under strict observation of this manual.
- This product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.

In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief description of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

**mindray**, **MINDRAY**, **迈瑞 MINDRAY** and **M** are the registered trademarks or trademarks owned by Mindray in China.

© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

### Product Introduction

#### Intended Use

The 551B veterinary SpO<sub>2</sub> sensor includes sensor pads, cables, connector and two clips with different size.

The reusable 551B veterinary SpO<sub>2</sub> sensor is applicable to canine and feline animals for spot check or continuous noninvasive measurement of arterial oxygen saturation and pulse rate.

The 551B veterinary SpO<sub>2</sub> sensor must be used in connection with the monitoring device (including accessories) specified by our company that supports the veterinary SpO<sub>2</sub> measurement function. The user shall be responsible for ensuring the compatibility between this sensor and the monitoring device (including accessories) by either referring to the instructions for using this sensor or contacting us before use. These sensors must be used by trained clinical professionals. Otherwise, it may damage the accessory or equipment, or bring harm to the animal.

Select the corresponding SpO<sub>2</sub> extension cable based on the type of SpO<sub>2</sub> connector.

#### Contraindications

Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or burns. Inspect the sensor site every two hours and move the sensor if the skin quality changes. Change the application site every four hours. For animals with poor peripheral blood circulation or sensitive skin, inspect the sensor site more frequently.

#### Warnings

- Check the sensor and the cable before using. Do not use them for SpO<sub>2</sub> monitoring in case of damage, block, deterioration or visible contamination.
- Change the clip in case of cracking, breakage or loss of spring tension.
- For disposal of accessories, follow local regulations or your hospital's policy regarding disposal of such accessories. Do not dispose randomly.
- To reduce the impact on the measurement of electrical surgical equipment, do not place the sensor cable near to the power lines of surgical equipment.
- The SpO<sub>2</sub> sensor must be used in connection with the specified equipment only. The user shall be responsible for ensuring the compatibility between this sensor, the extension cable and other equipment by either referring to the instructions for using this sensor or contacting us before use.
- Carefully route cables to reduce the possibility of animal entanglement or strangulation.
- Do not bend the cable excessively for a prolonged time.
- Do not use SpO<sub>2</sub> sensor during MRI scanning. Conducted current may cause burns.
- If using SpO<sub>2</sub> sensor during full body irradiation, keep the sensor out of the irradiation field. If sensor is exposed to the irradiation, the reading might be inaccurate or the product might read zero for the duration of the active irradiation period.
- Do not place the product in an environment that exceeds the stated range.
- Do not modify the SpO<sub>2</sub> sensor in any way. Modification may affect performance and accuracy.

#### Notes

- Performance of the sensor may be compromised by motion. Try to keep the animal still and avoid the measured site suffering excessive motion.
- Always choose a site that is well perfused and will completely cover the sensor's detector window. If the sensor site has very low perfusion, the SpO<sub>2</sub> reading may be inaccurate. As tongue's perfusion is generally higher than the ears and other parts, so tongue is the preferred sensor site when performing the measurement.
- A proper sensor clip shall be selected according to different animals and application sites. Large and small clip can be adapted to different animals and application sites.
- Tight sensor may cause venous pulsation, obstructed blood circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurement, while loose sensor may lead to erroneous optical alignment or falling off. If the sensor is too tight because the application site is too large or becomes too large due to edema, excessive pressure may result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.
- Inspect the application site every two hours and change the application site when taking a long time continuous monitoring. Some animals may require more frequent inspections, such as animals with perfusion disorder or sensitive skin. Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering, burns, venous pulsation, obstructed blood circulation, pressure marks, pressure necrosis and so on.
- Using the sensor in the presence of bright lights may result in inaccurate measurements. In such case, cover the sensor with an opaque material.

#### Applying SpO<sub>2</sub> sensor

Be sure to understand all warnings listed in the previous section before applying any SpO<sub>2</sub> sensor to the animal. Before use, verify whether the SpO<sub>2</sub> sensor is intact. In case of any damage, do not apply it to the animal.

##### Applying SpO<sub>2</sub> sensor

1. Clean the sensor and clip separately before and after each use.
2. Open the clip by pressing with the thumb and forefinger. Place one of the sensor pad into the clip slot. Adjust the sensor pad until the sensor pad is fully engaged in the clip slot, as shown in Figure 1. So it is with another sensor pad. The sensor clip is now ready to be applied to the animals, as shown in Figure 2.  
**Note: Be sure that the sensor pads are oriented as shown in the illustrations. The optical components must face each other directly.**
3. Place the clip to the measurement site. Ensure that the sensor's optical components positioned on the center of the tongue, as shown in Figure 3.  
**Note: The preferred application site for canine and feline animals is tongue. The other sites such as ear, lip, toe, prepuce or vulva can also be measured.**
4. Make sure that the sensor cable is positioned along the side of the animal's face and body to avoid entanglement with the animal, as shown in Figure 4.
5. Connect the extension cable and the sensor. Then connect the extension cable to the pulse oximetry equipment by referring to the operator's manual of the equipment.


Figure 1


Figure 2

Figure 3

Figure 4

#### Checking the Functionality of the Sensor

Replace the sensor if the following situations occur when the sensor is connected to the instrument:

- The "Sensor Fault" alarm occurs.
- There is no "Sensor Off" alarm when the sensor is closed but not applied to the ear lobe.
- The "Sensor Off" alarm occurs when the sensor is applied to the ear lobe.

Start SpO<sub>2</sub> monitoring only when you are sure that the sensor functions normally.

#### Removing SpO<sub>2</sub> Sensor

For the comfort of the animal and to avoid damaging the sensor do not pull on the cable when removing the sensor and clip.

1. Open the clip and remove the sensor and clip from the animal.
2. Hold and release the sensor pad from the clip slot.

#### Notes

- Failure to apply the sensor properly may cause inaccurate measurements. If the sensor cannot detect the pulse signal, it may be incorrectly positioned, or the application site is too thick or thin. If any of the situations occurs, reposition the sensor or try another sensor site. If the application site is covered with fur, shave the site and reapply the sensor.
- Do not use supplemental tape to adhere the clip and sensor directly to the application site, this can restrict blood flow and cause inaccurate measurements.
- To verify correct placement of the sensor on the site, manually take a pulse rate reading from another site on the animal if the sensor is placed correctly, the manual measurement should correspond to the pulse rate reading displayed on the monitor.

#### Measurement Limitations

If you doubt the measured SpO<sub>2</sub>, check animal's vital signs first. Then check the monitor and SpO<sub>2</sub> sensor.

The following factors may influence the accuracy of measurement:

- Ambient light
- Physical movement (animal and imposed motion)
- Diagnostic testing
- Low perfusion
- Electromagnetic interference, such as MRI environment
- Electrosurgical units
- Dysfunctional haemoglobin, such as carboxyhemoglobin (COHb) and methemoglobin (MetHb)
- Presence of certain dyes, such as methylene and indigo carmine
- Inappropriate positioning of the SpO<sub>2</sub> sensor, or use of incorrect SpO<sub>2</sub> sensor
- Drop of arterial blood flow to immeasurable level caused by shock, anemia, low temperature or vasoconstrictor.

#### Cleaning/Disinfection

##### Cautions

- Do not immerse the connector of the sensor in water or any of the cleaning solutions or disinfectants.
- Exercise caution during cleaning/disinfection to avoid wetting the pins.
- We recommend that the sensor be disinfected only when necessary as determined by your hospital's policy, to avoid long term damage to the sensor.
- Use cleaning solutions or disinfectants recommended by our company. Clean or disinfect the sensor according to the method described in operator's manual, otherwise, it might affect product performance or life.

##### Cleaning

1. Clean the sensor or clip with cotton or soft cloth moistened with water.
2. After cleaning, wipe off the water with a soft cloth.
3. Allow the sensor or clip to air dry.

##### Disinfection

The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde-type (2%) liquid disinfectants (using glutaraldehyde-type (2%) liquid disinfectants to disinfect the clip, corrosion inhibitor must be added and blending).

1. Clean the sensor or clip as instructed above.
2. Disinfect the sensor or clip with cotton or soft cloth moistened with one of the recommended disinfectants.
3. After disinfection, be sure to wipe off the disinfectant left on the sensor or clip with a soft cloth moistened with water.
4. Allow the sensor or clip to air dry.

#### Specifications

Cable Length	1.3m
Wavelength Range	600nm - 1000nm
Maximum photic output	18 mW
Note: The information about the wavelength range and maximum photic output consumption can be especially useful to clinicians (for example, when photodynamic therapy is performed).	

	SpO <sub>2</sub>	Pulse rate
Measurement range	0% - 100%	20 - 254 bpm
Accuracy	70 to 100%: ±3% 0% to 69%: Not specified.	±3 bpm
Environmental specifications	Temperature	Relative humidity
Operating conditions	0°C - 50°C	15% - 95% (noncondensing)
Storage conditions	-30°C - 70°C	10% - 95% (noncondensing)
		Barometric
		57.0kPa - 107.4kPa
		16.0kPa - 107.4kPa

#### Symbols

	Refer to instruction manual/booklet		Latex free
	Date of manufacture		Manufacturer
	Batch code		Part number
	Comply with the requirements of the Council Directive 93/42/EEC (Medical Device Directive).		Protection against vertically falling water drops when ENCLOSURE tilted up to 15°
	Authorized representative in the European Community		Federal law (USA) restricts this device to sale by or on the order of a physician.
	Dispose of in accordance to your country's requirements		

#### Company Contact

Website: [www.mindray.com](http://www.mindray.com)

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen 518057 P.R. China

Tel: +86 755 81888998

Fax: +86 755 26582680

E-mail: [service@mindray.com](mailto:service@mindray.com)

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, 20537 Hamburg, Germany

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