


R-LOCK (FLL, MLL)



1. Adapter with integrated CODAN Check Valve, with Male and Female Luer-Lock connector, prevents backflow of solution. Protection cap and the plug maintain the sterility of the set.
 - Minimum closure flow: 0,1 ml/h
 - Flow reduction: 34%
 - Opening pressure: 5 cm water column
 - Priming: without airbubbles
 - Pressure tightness: according to ISO 8536 part 12
 - Residual volume: 0,15 ml.

Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times. Subject to alterations without prior notice!

| REF | Product Description | Packaging Units |
|---------|---------------------|-----------------|
| 16.5250 | R-LOCK | 150/2400/57600 |