

325

AGA

Minigas







ARTEMA

REF 60-12001-00 Calibration Gas

5% CO₂ ±0.03%_{abs}; 1.5% ISO ±0.15%_{abs}

45% O₂ ±0.5%_{rel}; balance N₂O (nom. 48.5%)

 Only for test. Not for clinical use.



AGA

Minigas

WARNI

Test Gas

Made by AGA Gas AB for

Shenzhen Mindray Bio-Medical Electronics Co. Ltd

1,5% ISOFLURANE ±0.15% abs
5% CARBON DIOXIDE ±0.03% abs
45% OXYGEN ±0.5% rel
Balance DINITROGEN MONOXIDE

Filling date: 2015-09-24

Material number: 112462

Production order: 100392161

Certificate number: 100375139

Expiry date: 2016-03-24

Risk phrase:

May cause or intensify fire; oxidizer.

Safety phrase:

Keep away from open flames – No smoking.

Pressurized container: Do not pierce or burn, even after use.

Pressurized container for one time use.

Do not expose to temperatures exceeding +50°C/122

UN 1950 AEROSOL

Return emtied container to AGA Gas AB for recyclin

NOT REFILLABLE